




Safety ECOLOGICAL Funnel

Does your open waste container make you Fume?



Complete Systems



Safe
Effective
Economical



"Simple yet efficient"

Safety ECO Funnel Purpose & Benefits

- Safety ECO Funnel solves "the open waste container problem".
- Safety ECO Funnel prevents emission of 99.9% of volatile waste solvents, acids, and bases in the laboratory.
- Safety ECO Funnel is a cost effective alternative to EPA, OSHA and Fire Department fines and penalties.
- Significantly reduces chemical evaporation and environmental pollution of toxic and hazardous compounds in the laboratory.
- Unique patented product meets most local, state and federal safety requirements and regulations.
- Improves laboratory safety and reduces fire hazard.
- Safety ECO Funnel is reusable and has an average cost of approximately £0.04 per day assuming a four-year product use.
- Safety ECO Funnel protects laboratory scientists from exposure to toxic "fat soluble" compounds such as Dichloromethane, Hexane, Heptane, Acetonitrile, Disulfides etc.
- Safety ECO Funnel is recognized as a state of the art product in laboratory safety by environmental health and safety officers and industry professionals.
- Safety ECO Funnel is widely used by Pharmaceutical, Biotech, Academic and Government laboratories.

ORDERING INFORMATION



CalPacLab Catalog

EF-4-38C-System

Description

Consist of a Safety **ECO** Funnel with cover, LDPE bottle & Secondary Container to hold entire contents of funnel and waste container in case of a spill (see components listed below)

Price

£93.00

*Components | EF-4-38C
EC-2004
SC-3004CPB

4" Safety **ECO** Funnel, 38-430 screw cap, HDPE
4L Narrow-Mouth Bottle, 38-430 screw cap, LDPE
4L Secondary Container with the base, HDPE

£52.00

£29.00

£23.00



EF-3004C-System

Consist of a Safety **ECO** Funnel with cover, LDPE bottle & Secondary Container to hold entire contents of funnel and waste container in case of a spill (see components listed below)

£93.00

*Components | EF-3004C
EC-2004
SC-3004CPB

8" Safety **ECO** Funnel, 38-430 screw cap, HDPE
4L Narrow-Mouth Bottle, 38-430 screw cap, LDPE
4L Secondary Container with the base, HDPE

£52.00

£29.00

£23.00



EF-8-Justrite-B2-System

Justrite® Safety **ECO** Funnel & 2 Gal translucent Justrite® Container (see components listed below)

£215.00

EF-8-Justrite-B5-System

Justrite® Safety **ECO** Funnel & 5 Gal translucent Justrite® Container (see components listed below)

£189.00

*Components | EF-8-Justrite-B
EC-12751
EC-12754

8" Safety **ECO** Funnel, HDPE
2 Gal translucent, FM Approved, Fire rated safety can, polyethylene base can
5 Gal translucent, FM Approved, Fire rated safety can, polyethylene base can

£92.00

£137.00

£137.00

*Components also sold separately

Safety **ECO** Funnel can also be customized to your desired container and screw cap size.

EF-4-38C-System

- Top diameter: 4" with positive latch
- Screw Cap: 38-430
- Stem Length: 11" long
- Funnel Capacity: 250 mL
- Recommended Bottle: 4 L fits any 38-430 screw cap Safety **ECO** Funnel
- Complete System sold with EC-2004 white LDPE bottle and SC-3004CPB HDPE Secondary Container

EF-3004C-System

- Top diameter: 8" with positive latch
- Screw Cap: 38-430
- Stem Length: 11" long
- Funnel Capacity: 1750 mL
- Recommended Bottle: 4 L fits any 38-430 screw cap Safety **ECO** Funnel
- Complete System sold with EC-2004 white LDPE bottle and SC-3004CPB HDPE Secondary Container

EF-8-Justrite-B2-System & EF-8-Justrite-B5-System

- Top diameter: 8" with positive latch
- Funnel Capacity: 1750 mL
- Recommended Container: 2 & 5 Gal translucent Justrite Can®

Please visit our website at www.presearch.co.uk

Safety **ECO** Funnel Technical Specifications

- Positive latch
- Gasket in the lid
- Industrial strength screw cap
- Chemically resistant HDPE funnel
- Filter to catch large solid particles or stir bars
- Safety **ECO** Funnel are also available with an HPLC hook up fitting (*optional*)



Call to ask our salesperson if your company qualifies for a free demonstration and display Safety **ECO** Funnel System.

UK Distributor: Presearch Ltd
System House, 59-61 Knowl Piece, Hitchin
Herts SG4 0TY
T: 01462 442424 www.presearch.co.uk
F: 01462 442429 enquiries@presearch.co.uk



EVOLUTION OF FUNNELS
WWW.CALPACLAB.COM
1-888-322-5722

Prices are current and may change without notice.
Note: Availability for certain products may be limited by federal, state, provincial, or local licensing requirements. All prices printed are in GBP.

Justrite® is a registered trademark of Justrite Mfg. Co., L.L.C.